



Innovative Cryo Enterprises L.L.C.

I.C.E. Oocyte Vitrification Kit

Vitrification Media V1, V2, V3



Innovative Cryo Enterprises L.L.C.

ICE Oocyte Vitrification Instructions For Use Testing and Cautions

Innovative Cryo Enterprises LLC
317 Springfield Road
Linden, New Jersey 07036
USA
973-632-8635

Intended Use:

To be used for vitrification of human oocytes.
Not to be used with blastocysts.

Quality Control Testing:

Sterility (SAL 10^{-3})

pH

Endotoxin Tested ≤ 0.5 EU/ml (USP)

Mouse Embryo assay (MEA)

Note: The results of each batch are listed on a Certificate of Analysis, which is available upon request.

Storage instructions and stability:

Store in original container at 2°C to 8°C protected from light.

The product is packaged in tubes and can be used multiple times.

Keep media at 2°C to 8°C when not in use, do not keep at room temperature for extended periods of time.

Precautions and Warnings:

Do not use product if:

- 1) Product packaging appears damaged or if the seal is broken.
- 2) Expiration date has been exceeded.
- 3) Any solution becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

To avoid contamination problems, handle using aseptic techniques.

Keep away from sunlight.

Consult operating instructions.

ICE Oocyte vitrification and thaw solutions contain the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure the patient is not sensitized to this antibiotic.

The long term safety of oocyte vitrification on children born using this technique is unknown.

This product contains albumin, a product of human blood.

*Human source materials used in the manufacture of this product have been tested with FDA licensed kits, and found to be non-reactive to the antibodies for Hepatitis B surface antigen (HsbAg), antibodies to Hepatitis C (HCV) and antibodies to Human Immuno-deficiency Virus (HIV). Donors of the source material have also been screened for CJD. However, no test method offers complete assurance that

jstachecki@gmail.com 973-632-8635 www.icevitrification.com



Innovative Cryo Enterprises L.L.C.

products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal precautions.

Caution: US federal law restricts this device to sale by or on the order of a physician (Rx only) or a practitioner trained in its use.

Caution: The user should read and understand the Directions for Use, Warnings, and Precautions, and be trained in the correct procedure before using the ICE Oocyte Vitrification kits for storage of human oocytes.

Quality Assurance

All solutions are membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10^{-3} .

Each lot of ICE Oocyte Vitrification & Thaw receives the following tests:

pH

Endotoxin by LAL methodology

Biocompatibility and functionality by mouse embryo assay (MEA)

Sterility by the current USP Sterility Test <71>



I.C.E. OOCYTE VITRIFICATION PROTOCOL

PART 1

Materials Required:

- ICE Oocyte Vitrification Media V1, V2, & V3

Materials required but not included:

- 0.25cc or 0.5cc sterile cryo-straw(s) or micro-volume device (cryo-top, cryo-tec, cryo-loc, etc.)
- Heat sealer
- Container for liquid nitrogen
- Storage tank
- Sterile petri dishes (60x9mm or 4-6 well dish or similar)
- 1ml sterile syringe
- Stopwatch or Timer
- Liquid nitrogen
- Micropipettes for moving cells (270-300um inner diameter)
- Storage cane and goblet
- Stereomicroscope (Heating plate off)

CAUTION: Use for vitrification of oocytes only. Do not freeze blastocysts.

I.C.E. OOCYTE VITRIFICATION PROTOCOL

PART 2

Preparation for Vitrification:

- 1) Bring V1, V2, & V3 to room temperature (23-27°C).
- 2) Label Cryostraw (0.25cc) or storage device with necessary patient information. (1-5 oocytes per storage device is recommended).
- 3) Fill container with liquid nitrogen.
- 4) Turn on heat sealer and test. (If using straws).
- 5) Label one petri dish or 4-6 well plate: V1, V2, V3 for each straw or device to be stored. (See Figure 2.1).






Figure 2.1

I.C.E. OOCYTE VITRIFICATION PROTOCOL

Part 3. Oocyte Vitrification


Note: All steps are performed at 23-26°C. Use a timer for all steps.

- 1) Pipet 75-100ul of V1, V2, & V3 into the labeled dish and replace the cover. (Do not use oil!)
- 2) Collect and transfer 1-5 oocytes with minimal media to the V1 drop. *The number of oocytes will depend on how many you will vitrify per device.*
- 3) Incubate in **V1** for **5 min.**
- 4) Transfer the oocytes with minimal media to **V2** for **2.5 min.**
- 5) Transfer the oocytes with minimal media to **V3** for **1-2 min.** Move oocytes to a clean area of the drop several times to ensure complete mixing. Move the oocytes to a small drop of 10-15ul, and load the entire drop into the straw. *Be sure to soak the oocytes at least 2 min in V3 prior to storing in LN2.*




Instructions for loading and storing oocytes in 0.25cc or 0.5cc straws.

6) Load oocytes in a 0.25cc or 0.5cc straw between columns of V3. See Figure 3.3 & 3.4. *Aspirate the first column of V3 to the cotton plug, as this will allow the columns to remain intact.*



7) Heat seal both ends of the straw. See Figure 3.4. *You must heat-seal both ends of the straw! Periodically check seal under microscope to ensure a good seal.*




8) Plunge and store in LN₂.

Instructions for loading and storing oocytes on a micro-volume device.

6) Load oocytes onto micro-volume device (cryo-top, cryo-tec, rapid-i, cryo-loc, etc.). Be careful not to load too much media onto device. See Figure 3.5.

Always refer to the manufacturers recommendations for loading and storing the device.



7) Plunge and store in LN₂.

For optimal performance:

Wait 1 hr after egg collection, strip most of the cumulus cells off, wait 1 hr then vitrify.

Leaving some cumulus cell on helps with locating the eggs upon thawing.

Label the straw directly or use a flattened 0.5cc straw as a label that can be heat-sealed onto the 0.25cc straw.

Push the cotton plug in 2cm to have room for heat-sealing.

Use a 60 mm dish to freeze in. Use large drops (75-100 µl) with no oil overlay.

Because the time allowed in V3 can be up to 5 min, there is time to prepare multiple straws for vitrification.

Use a new dish (new medium) for each set of cells frozen.

Aspirate the first column of V3 to the cotton plug, as this will allow the columns to remain intact.

Store media at 2-8°C.



Innovative Cryo Enterprises L.L.C.

I.C.E. Oocyte Thaw Kit

Thaw Media T1, T2, T3, T4, T5



Innovative Cryo Enterprises L.L.C.

ICE Oocyte Thaw Instructions For Use Testing and Cautions

Innovative Cryo Enterprises LLC
317 Springfield Road
Linden, New Jersey 07036
USA
973-632-8635

Intended Use:

To be used for thawing oocytes frozen using ICE Oocyte Vitrification.
Not to be used with blastocysts.

Quality Control Testing:

Sterility (SAL 10^{-3})

pH

Endotoxin Tested ≤ 0.5 EU/ml (USP)

Mouse Embryo assay (MEA)

Note: The results of each batch are listed on a Certificate of Analysis, which is available upon request.

Storage instructions and stability:

Store in original container at 2°C to 8°C protected from light.

The product is packaged in tubes and can be used multiple times.

Keep media at 2°C to 8°C when not in use, do not keep at room temperature for extended periods of time.

Precautions and Warnings:

Do not use product if:

- 1) Product packaging appears damaged or if the seal is broken.
- 2) Expiration date has been exceeded.
- 3) Any solution becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

To avoid contamination problems, handle using aseptic techniques.

Keep away from sunlight.

Consult operating instructions.

ICE Oocyte Thaw solutions contain the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure the patient is not sensitized to this antibiotic.

The long term safety of oocyte vitrification on children born using this technique is unknown.

This product contains albumin, a product of human blood.

*Human source materials used in the manufacture of this product have been tested with FDA licensed kits, and found to be non-reactive to the antibodies for Hepatitis B surface antigen (HsbAg), antibodies to Hepatitis C (HCV) and antibodies to Human Immuno-deficiency Virus (HIV). Donors of the source material have also been screened for CJD. However, no test method offers complete assurance that

jstachecki@gmail.com 973-632-8635 www.icevitrification.com



Innovative Cryo Enterprises L.L.C.

products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal precautions.

Caution: US federal law restricts this device to sale by or on the order of a physician (Rx only) or a practitioner trained in its use.

Caution: The user should read and understand the Directions for Use, Warnings, and Precautions, and be trained in the correct procedure before using the ICE Oocyte Thaw kits for recovery of vitrified human oocytes.

Quality Assurance

All solutions are membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10^{-3} .

Each lot of ICE Oocyte Vitrification & Thaw receives the following tests:

pH

Endotoxin by LAL methodology

Biocompatibility and functionality by mouse embryo assay (MEA)

Sterility by the current USP Sterility Test <71>



I.C.E. OOCYTE THAW PROTOCOL

PART 1

Materials Required:

- ICE Oocyte Thaw Media T1, T2, T3, T4, & T5

Materials required but not included:

- Container for liquid nitrogen
- Sterile petri dishes (60x9mm or 4-6 well dish or similar)
- 1ml sterile syringe
- Stopwatch or Timer
- Liquid nitrogen
- Micropipettes for moving cells (270-300um inner diameter)
- Stereomicroscope (Heating plate off)
- Disposable gloves
- Scissors (sterile)
- 30°C water bath
- Mineral oil
- Culture medium appropriate for oocyte stage cells

I.C.E. OOCYTE THAW PROTOCOL

PART 2

Preparation for Thawing:

1) The day before thawing prepare dishes with culture medium for the thawed oocytes and allow to equilibrate in incubator overnight. It is recommended to culture thawed cells in media containing 20% HSA or a similar protein, after thawing and before ICSI and after ICSI until fertilization check.

2) Label one petri dish or 6 well plate: T1, T2, T3, T4, T5. Do this for each straw or device to be thawed. (See Figure 2.1).

3) Pipet 75-100ul of T1, T2, T3, T4, T5 into the labeled dish and cover the drops with mineral oil. (Return media to refrigerator at 2-8°C).

4) Fill container with liquid nitrogen.

5) Select straw(s) or devices to be thawed and place into container with liquid nitrogen.

6) Bring T1 to room temperature (23-27°C).






Figure 2.1

I.C.E. OOCYTE THAW PROTOCOL

PART 3 Oocyte Thawing

Please Note: All steps are performed at 23-26°C unless otherwise indicated. Use a timer for all steps. Oocytes are most often vitrified in a conventional 0.25cc straw or 0.5cc straw, but other micro-volume devices can be used. The cells are located in the middle column of a straw or the tip of the micro-volume device. See Figure 3.1.





For thawing straws:

1) Label one petri dish T1 and pipet 250-300ul T1 into dish. Do not cover with oil. See Figure 3.2.

2) Remove straw from liquid nitrogen and hold for **3 sec in room temperature air**. Use a timer! See Figure 3.4.


3) Submerge straw completely in **37°C water for 7 sec**. Figure 3.4.

4) Wipe excess water off straw and cut off heat-sealed end of the straw. Insert syringe and cut off the other heat-sealed end. See Figure 3.5.



5) Empty contents of straw into a 60mm petri dish. See Figure 3.6.


6) Watch through the microscope to see the oocytes come out of the straw. They will appear flat and "ghost" like. You may see what is an empty zona or dead oocyte. The oocyte is most likely alive despite initial appearance.



Go to Step #7.

For thawing micro-volume devices (cryo-top, cryo-loc, cryo-tec, etc):

1) Label one small petri or organ culture dish T1 and pipet 500ul-1500ul T1 into dish. Do not cover with oil. Place lid on dish. Incubate dish at 37-38°C for 10-15 min to warm the media. Make sure T1 media is warm and at 37-38°C before using.




2) Remove micro-volume device from liquid nitrogen and immediately submerge the end containing the oocytes completely into 37-38°C T1. Wait until the oocytes fall off the device. See Figure 3.7. Always refer to micro-volume device manufacturers recommendations for thawing.

3) Go to Step #7.

All steps are performed at 23-26°C in microdrops under oil.

- 7) Transfer the oocytes to T1 for 2 min. See Figure 3.8.
- 8) Transfer the oocytes to T2 for 3 min.
- 9) Transfer the oocytes to T3 for 3 min.
- 10) Transfer the oocytes to T4 for 3 min.
- 11) Transfer the oocytes to T5 and place dish on a 37°C warmer for 3 min.
- 12) Transfer the oocytes to equilibrated and warmed culture medium in a 37°C incubator. Culture 2-3h before ICSI.



For optimal performance:

For best results use a new dish (new medium) for each straw thawed.
 Culture thawed cells in media containing 20% HSA for the first 24 hours after thawing.
 After the first 24 hours, you can switch to your normal protein concentration for culturing (no less than 10% HSA is recommended).
 Store media at 2-8°C.